

AIR-OXYGEN % CHART: Choose Desired % - Choose Desired Total Flow - Set Air and Oxygen LPM Accordingly

O2 %	T-FLOW	LPM O2	LPM AIR
25%	9	0.5	8.5
	9.5	0.5	9
	10	0.5	9.5
	10.5	0.5	10
26%	7.5	0.5	7
	8	0.5	7.5
	8.5	0.5	8
27%	6.5	0.5	6
	7	0.5	6.5
28%	5.5	0.5	5
	6	0.5	5.5
	11	1	10
29%	5	0.5	4.5
	9.5	1	8.5
	10	1	9
	10.5	1	9.5
30%	4.5	0.5	4
	8.5	1	7.5
	9	1	8
31%	4	0.5	3.5
	8	1	7
	11.5	1.5	10
32%	3.5	0.5	3
	7	1	6
	7.5	1	6.5
	10.5	1.5	9
	11	1.5	9.5
33%	6.5	1	5.5
	9.5	1.5	8
	10	1.5	8.5
34%	3	0.5	2.5
	6	1	5
	9	1.5	7.5
35%	12	2	10
	5.5	1	4.5
	8.5	1.5	7
	11	2	9
36%	11.5	2	9.5
	8	1.5	6.5
	10.5	2	8.5
37%	2.5	0.5	2
	5	1	4
	7.5	1.5	6
	10	2	8
	12	2.5	9.5
38%	12.5	2.5	10
	7	1.5	5.5
	9.5	2	7.5
39%	11.5	2.5	9
	4.5	1	3.5
	6.5	1.5	5
	9	2	7
	11	2.5	8.5
40%	13	3	10
	8.5	2	6.5
	10.5	2.5	8
41%	12.5	3	9.5
	4	1	3
	6	1.5	4.5
	8	2	6
	10	2.5	7.5
	12	3	9
13.5	3.5	10	
	7.5	2	5.5

O2 %	T-FLOW	LPM O2	LPM AIR
42%	9.5	2.5	7
	11.5	3	8.5
	13	3.5	9.5
43%	5.5	1.5	4
	9	2.5	6.5
	11	3	8
	12.5	3.5	9
44%	3.5	1	2.5
	7	2	5
	8.5	2.5	6
	10.5	3	7.5
	12	3.5	8.5
45%	13.5	4	9.5
	14	4	10
	5	1.5	3.5
	6.5	2	4.5
46%	10	3	7
	11.5	3.5	8
	13	4	9
	8	2.5	5.5
47%	9.5	3	6.5
	11	3.5	7.5
	12.5	4	8.5
	14	4.5	9.5
	14.5	4.5	10
48%	3	1	2
	4.5	1.5	3
	6	2	4
	7.5	2.5	5
	9	3	6
	10.5	3.5	7
	12	4	8
	13.5	4.5	9
15	5	10	
49%	11.5	4	7.5
	13	4.5	8.5
	14.5	5	9.5
	7	2.5	4.5
50%	8.5	3	5.5
	10	3.5	6.5
	12.5	4.5	8
	14	5	9
	15.5	5.5	10
51%	5.5	2	3.5
	9.5	3.5	6
	11	4	7
	13.5	5	8.5
52%	15	5.5	9.5
	4	1.5	2.5
	6.5	2.5	4
	8	3	5
	10.5	4	6.5
	12	4.5	7.5
	13	5	8
53%	14.5	5.5	9
	16	6	10
	9	3.5	5.5
	11.5	4.5	7
	14	5.5	8.5
	16.5	6.5	10
54%	2.5	1	1.5
	5	2	3
	7.5	3	4.5
	10	4	6
	12.5	5	7.5
55%	15	6	9
	17.5	7.5	10
	3.5	1.5	2
	7	3	4
56%	10.5	4.5	6.5
	11.5	5	6.5
	14	6	8
	15	6.5	8.5
	16.5	7	9.5
	17.5	7.5	10
	4.5	2	2.5
	8	3.5	4.5
57%	9	4	5
	12.5	5.5	7
	13.5	6	7.5
	14.5	6.5	8
	16	7	9
	17	7.5	9.5
	18	8	10
	5.5	2.5	3
58%	6.5	3	3.5
	10	4.5	5.5
	11	5	6
	12	5.5	6.5
	13	6	7
	15.5	7	8.5
	16.5	7.5	9
59%	17.5	8	9.5
	7.5	3.5	4
	8.5	4	4.5
	9.5	4.5	5
	14	6.5	7.5
	15	7	8
	16	7.5	8.5
60%	17	8	9
	18	8.5	9.5
	10.5	5	5.5
	11.5	5.5	6
	12.5	6	6.5
	13.5	6.5	7
	14.5	7	7.5
61%	15.5	7.5	8
	16.5	8	8.5
	17.5	8.5	9
	6	3	3
	8	4	4
	9	4.5	4.5
	10	5	5
62%	11	5.5	5.5
	12	6	6
	13	6.5	6.5
	15	7.5	7.5

AIR-OXYGEN % CHART: Choose Desired % - Choose Desired Total Flow - Set Air and Oxygen LPM Accordingly

O2 %	T-FLOW	LPM O2	LPM AIR
61%	3	1.5	1.5
	4	2	2
	5	2.5	2.5
	7	3.5	3.5
	14	7	7
	16	8	8
	17	8.5	8.5
18	9	9	
62%	10.5	5.5	5
	11.5	6	5.5
	12.5	6.5	6
	13.5	7	6.5
	14.5	7.5	7
	15.5	8	7.5
	16.5	8.5	8
17.5	9	8.5	
63%	7.5	4	3.5
	8.5	4.5	4
	9.5	5	4.5
	14	7.5	6.5
	15	8	7
	16	8.5	7.5
	17	9	8
18	9.5	8.5	
64%	5.5	3	2.5
	6.5	3.5	3
	10	5.5	4.5
	11	6	5
	12	6.5	5.5
	13	7	6
	15.5	8.5	7
16.5	9	7.5	
17.5	9.5	8	
65%	4.5	2.5	2
	8	4.5	3.5
	9	5	4
	12.5	7	5.5
	13.5	7.5	6
	14.5	8	6.5
	16	9	7
17	9.5	7.5	
18	10	8	
66%	3.5	2	1.5
	7	4	3
	10.5	6	4.5
	11.5	6.5	5
	14	8	6
	15	8.5	6.5
	16.5	9.5	7
17.5	10	7.5	
67%	6	3.5	2.5
	8.5	5	3.5
	9.5	5.5	4
	12	7	5
	13	7.5	5.5
	14.5	8.5	6
	15.5	9	6.5
17	10	7	
68%	2.5	1.5	1
	5	3	2
	7.5	4.5	3
	10	6	4
	11	6.5	4.5
12.5	7.5	5	

O2 %	T-FLOW	LPM O2	LPM AIR
68%	13.5	8	5.5
	15	9	6
	16	9.5	6.5
69%	9	5.5	3.5
	11.5	7	4.5
	14	8.5	5.5
	15.5	9.5	6
	16.5	10	6.5
70%	4	2.5	1.5
	6.5	4	2.5
	8	5	3
	10.5	6.5	4
	12	7.5	4.5
	13	8	5
	14.5	9	5.5
16	10	6	
71%	5.5	3.5	2
	9.5	6	3.5
	11	7	4
	13.5	8.5	5
	15	9.5	5.5
72%	7	4.5	2.5
	8.5	5.5	3
	10	6.5	3.5
	12.5	8	4.5
	14	9	5
15.5	10	5.5	
73%	11.5	7.5	4
	13	8.5	4.5
	14.5	9.5	5
74%	3	2	1
	4.5	3	1.5
	6	4	2
	7.5	5	2.5
	9	6	3
	10.5	7	3.5
	12	8	4
13.5	9	4.5	
15	10	5	
75%	8	5.5	2.5
	9.5	6.5	3
	11	7.5	3.5
	12.5	8.5	4
	14	9.5	4.5
14.5	10	4.5	
76%	5	3.5	1.5
	6.5	4.5	2
	10	7	3
	11.5	8	3.5
	13	9	4
77%	3.5	2.5	1
	7	5	2
	8.5	6	2.5
	10.5	7.5	3
	12	8.5	3.5
	13.5	9.5	4
	14	10	4
78%	5.5	4	1.5
	9	6.5	2.5
	11	8	3
	12.5	9	3.5
79%	7.5	5.5	2
	9.5	7	2.5
	13	9.5	3.5

O2 %	T-FLOW	LPM O2	LPM AIR
80%	11.5	8.5	3
	4	3	1
	6	4.5	1.5
	8	6	2
	10	7.5	2.5
81%	12	9	3
	13.5	10	3.5
	8.5	6.5	2
82%	10.5	8	2.5
	12.5	9.5	3
	4.5	3.5	1
	6.5	5	1.5
	9	7	2
83%	11	8.5	2.5
	13	10	3
	7	5.5	1.5
84%	9.5	7.5	2
	11.5	9	2.5
	2.5	2	0.5
	5	4	1
	7.5	6	1.5
85%	10	8	2
	12	9.5	2.5
	12.5	10	2.5
86%	8	6.5	1.5
	10.5	8.5	2
	5.5	4.5	1
87%	8.5	7	1.5
	11	9	2
	11.5	9.5	2
88%	3	2.5	0.5
	6	5	1
	9	7.5	1.5
	12	10	2
89%	6.5	5.5	1
	9.5	8	1.5
	10	8.5	1.5
90%	3.5	3	0.5
	7	6	1
	7.5	6.5	1
	10.5	9	1.5
91%	11	9.5	1.5
	4	3.5	0.5
	8	7	1
92%	11.5	10	1.5
	4.5	4	0.5
	8.5	7.5	1
93%	9	8	1
	5	4.5	0.5
	9.5	8.5	1
94%	10	9	1
	10.5	9.5	1
	5.5	5	0.5
95%	6	5.5	0.5
	11	10	1
	6.5	6	0.5
96%	7	6.5	0.5
	7.5	7	0.5
	8	7.5	0.5
97%	8.5	8	0.5
	9	8.5	0.5
	9.5	9	0.5
	10	9.5	0.5
98%	10.5	10	0.5
	11	10.5	1