

DELUXE MRI E-CYLINDER REGULATORS



Ideal for Transport

➤ **Unique DISS Input Fitting**

Allows switching between gas sources without disrupting the gas supply to the ventilator. Allows tanks to be changed out while connected to the wall source (in the event the tank is depleted during transport).

➤ **Built-in Pressure Gauge**

Not exposed to possible damage.

➤ **Built-in 0-15 LPM Flow Meter (O₂ regulator only)**

Provides emergent access to oxygen flow if needed for a manual resuscitator bag.

➤ **Built-in Permanent O-ring Tank Washer**

No need to search for lost washers.

➤ **MRI Compatible**

➤ **Preset to 50 PSI (344.7 kPa)**

Ordering Information:

Deluxe O₂ MRI Regulator- Catalog #1005M

Deluxe Air MRI Regulator- Catalog #1006M