


IC-2A PEDIATRIC THROUGH ADULT MRI VENTILATOR



Portable, Dependable, MR Conditional

Bio-Med Devices' IC-2A pediatric through adult MRI transport ventilator has been an industry favorite for many years. It is gas operated with no electronic or magnetic parts so it can be placed directly in the MRI room.

- Transport
- Emergency Room
- Calibrated Controls
- Built-in PEEP
- SIMV
- CPAP

IC-2A

Pediatric through Adult MRI Ventilator

Specifications

Parameter	Range
Gas Supply	50+/-5 PSI (345 +/-35 kPa)
Inspiratory Time (calibrated)	0.4 to 2 seconds
Expiratory Time (calibrated)	0.5 to 4 seconds
Rate	1 1/3 to 66 BPM
I/E	Infinitely Adjustable
Flowrate (calibrated)	0 to 75 LPM
Tidal Volume	0 to over 3000 ml
PEEP/CPAP	0 to 25 cmH2O
Pressure Trigger	-0.5 to -20 cmH2O

Ventilation Modes

CMV | ASSIST CONTROL | SIMV | PEEP | CPAP
Manual | Pressure Ventilation

Physical Characteristics

Size: 3-3/8" x 6" x 10-1/4" (8.57 x 15.24 x 26.04 cm)
Weight: 8-1/2 lbs. (3.86 kg.)

Ordering Information

800IAMRI IC-2A MRI Ventilator

Included Accessories

Two Disposable Patient Circuits
Standard Pole Mount Bracket
Oxygen Supply Hose
High Pressure Oxygen Wye
Supply Elbow Adapters
Test Lung

Optional Accessories

MRI Heavy Duty Stand with Tilt Bracket
MRI High/Low Flow Blender
MRI Cylinder Brackets
MRI E-Cylinder Deluxe Regulator- Oxygen
MRI E-Cylinder Deluxe Regulator- Air
MRI Aluminum E-Cylinder- Oxygen
MRI Aluminum E-Cylinder- Air

IC-2A Complete Package
includes Air/Oxygen Blender, Heavy Duty Roll Stand,
Cylinder Brackets, Hose Kit and E-Cylinder Deluxe Regulators
Order as Part #800IAMBC
(E-Cylinders sold separately)



Bio-Med Devices, Inc.
61 Soundview Road
Guilford, CT 06437
U.S.A.
Telephone: 800-224-6633 or 203-458-0202
Fax: 203-458-0440
Website: www.biomeddevices.com
E-mail: custserv@biomeddevices.com



Specifications subject to change without notice.

Distributed by

REV0208