

# BIO-MED DEVICES INC.

**NEW!**

## **DELUXE MRI E-CYLINDER REGULATORS**

**Ideal for Transport**



### ➤ **Unique DISS Input Fitting**

- ✓ Allows switching between gas sources without disrupting the gas supply to the ventilator.
- ✓ Allows tanks to be changed out while connected to the wall source (in the event the tank is depleted during transport).

### ➤ **Built-in Pressure Gauge**

- ✓ Not exposed to possible damage.

### ➤ **Built-in 0-15 LPM Flow Meter (O<sub>2</sub> regulator only)**

- ✓ Provides emergent access to oxygen flow if needed for a manual resuscitator bag.

### ➤ **Built-in Permanent O-ring Tank Washer**

- ✓ No need to search for lost washers.

### ➤ **MRI Compatible**

### ➤ **Preset to 50 PSI (344.7 kPa)**

#### Ordering Information:

Deluxe O<sub>2</sub> MRI Regulator- Catalog #1005M

Deluxe Air MRI Regulator- Catalog #1006M

Bio-Med Devices, Inc • 1445 Boston Post Road • Guilford, CT 06437  
800-224-6633 • Fax 203-458-0440  
[www.biomeddevices.com](http://www.biomeddevices.com)