




BIO-MED
DEVICES

Deluxe MRI E-Cylinder Regulators Ideal for Transport



Unique DISS Input Fitting

Allows switching between gas sources without disrupting the gas supply to the ventilator.
Allows tanks to be changed out while connected to the wall source (in the event the tank is depleted during transport).

Built-in Pressure Gauge

Not exposed to possible damage.

Built-in 0-15 LPM Flow Meter (O₂ regulator only)

Provides emergent access to oxygen flow if needed for a manual resuscitator bag.

Built-in Permanent O-ring Tank Washer

No need to search for lost washers.

MRI Compatible

Preset to 50 PSI (344.7 kPa)

Ordering Information:

1005M Deluxe O₂ MRI Regulator
1006M Deluxe Air MRI Regulator

Bio-Med Devices Inc.

61 Soundview Road

Guilford, CT 06437

USA

Telephone: 1-800-224-6633 or 203-458-0202

Fax: 203-458-0440

Website: www.biomeddevices.com

E-mail: custserv@biomeddevices.com

International@biomeddevices.com

DDOC012-091616