

Air/Oxygen Blenders 2001



Part Number: 2001

Flow Range: 2 to 120 LPM

Input Port Fittings: DISS (for NIST add Part # 2007 to your order)

Output Port Fittings: DISS

Output Port Flow Range:

Left – Plugged

Bottom – 15-120 LPM (no bleed)

2-108 LPM (with bleed)

Right – 2-108 LPM (with bleed)

Bleed for Flows <15 LPM: 12 LPM

Size: 3.6W x 5.4H x 4.6D inches

9.2W x 13.7H x 11.7D cm

Weight: 3.2 lbs./1.45 kg

Bio-Med Devices Inc.
61 Soundview Road
Guilford, CT 06437 USA
Telephone: 1-800-224-6633 or 203-458-0202
Fax: 203-458-0440
www.biomeddevices.com
E-mail: custserv@biomeddevices.com