


# Air/Oxygen Blenders 2002F70D



**Part Number:** 2002F70D

**Flow Range:** 15 to 120LPM

**Input Port Fittings:** DISS (for NIST add Part # 2007 to your order)

**Output Port Fittings:** DISS

**Output Port Flow Range:**

Left – 15-70 LPM Flowmeter (white O<sub>2</sub> label)

Bottom – 15-120 LPM

Right – Plugged

**Bleed:** None

**Size:** 4.9W x 7.0H x 4.6D inches  
12.5W x 17.8H x 11.7D cm

**Weight:** 3.5 lbs./1.6Kg