


# Air/Oxygen Blenders 2003P



**Part Number:** 2003P

**Flow Range:** 0 to 30 LPM

**Input Port Fittings:** DISS (for NIST add Part # 2007 to your order)

**Output Port Fittings:** DISS

**Output Port Flow Range:**

Left – 0-15 LPM Flowmeter (green O<sub>2</sub> label) Bottom  
– 3-30 LPM (no bleed)

0-30 LPM (with bleed)

Right – Knob to switch bleed on & off

**Bleed for Flows <3 LPM:** 3 LPM

**Size:** 5.6W x 6.7H x 4.6D inches  
14.2W x 17.0H x 11.7D cm

**Weight:** 3.5 lbs./1.6Kg