

Air/Oxygen Blenders 2004

Part Number: 2004

Flow Range: 2 to 120 LPM

Input Port Fittings: DISS (for NIST add Part # 2007 to your order)

Output Port Fittings: DISS

Output Port Flow Range:

Left – 15-120 LPM (no bleed)/ 2-108 LPM (with bleed)

Bottom – 15-120 LPM (no bleed)/ 2-108 LPM (with bleed)

Right – 2-108 LPM (with bleed)

Bleed for Flows <15 LPM: 12 LPM

Size: 4.8W x 5.4H x 4.6D inches

12.2W x 13.7H x 11.7D cm

Weight: 3.3 lbs./1.5 kg

