

Air/Oxygen Blenders 20047015EU

PartNumber:20047015EU

Flow Range: 2 to 120 LPM

Input Port Fittings: DISS (for NIST add Part # 2007 to your order)

Output Port Fittings: DISS

Output Port Flow Range:

Left – 15-70 LPM Flowmeter (green O2 label)

Bottom – 2-120LPM

0-30 LPM (with bleed)

Right – 1-15 LPM Flowmeter (green O2 label)

Bleed: 12 LPM

Size: 8.1W x 7.0H x 4.6D inches
20.6W x 17.8H x 11.7D cm

Weight: 4.1 lbs./ 1.86 kg

