

# Air/Oxygen Blenders

## 2004 F7015



**PartNumber:**2004F7015

**Flow Range:** 2 to 120 LPM

**Input Port Fittings:** DISS (for NIST add Part # 2007 to your order)

**Output Port Fittings:** DISS

**Output Port Flow Range:**

Left – 15-70 LPM Flowmeter (green O2 label)

Bottom – 15-120LPM

2-108 LPM (with bleed)

Right – 1-15 LPM Flowmeter (green O2 label)

**Bleed:** 12 LPM

**Size:** 8.1W x 7.0H x 4.6D inches  
20.6W x 17.8H x 11.7D cm

**Weight:** 4.1 lbs./ 1.86 kg